

# Instructions for wheelchair measurements

Arms of Aloha LLC

Need help? Call Lucy 808-435-3006 ext 6

Don't worry, we'll double-check the measurements for you, but this will get us in the ballpark so we can give you an estimated cost and bring the right equipment along.

**Fold of Flank:** with your dog lying on its side, measure from the toes to the fold of the flank: this is the flap of skin that connects the dog's upper leg to their body (like the webbing between your index finger and thumb). Leave your dog's leg in a natural position, slightly bent, instead of pulling it straight.

RESULT: \_\_\_\_\_



**Armpit to Rump:** Measure from the back of the front leg to the back of the rear leg. We only need this measurement for dogs 25 pounds or smaller.

RESULT: \_\_\_\_\_